

PRODUCT DATA

PRODUCT: DST/F

ISSUE NUMBER: 180707/3

<u>SUBSTANCE</u>		<u>THICKNESS</u>	
Adhesive	37gsm		
Tissue	12gsm	Adhesive & Carrier	85 micron
Adhesive	37gsm		
Liner	65gsm	Liner	55 micron

PRODUCT DESCRIPTION

A permanent tape with a tissue carrier. The adhesive is an aggressive rubber based system and the liner a two side release system which overlaps at each edge to give a 'fingerlift' for easy liner removal

APPLICATIONS

The original and industry standard fingerlift tissue tape, DST/F has been used for continuous stationery, mailing envelopes and padded envelopes, studio / commercial artwork, photographic, Point of Sale and splicing. Especially useful for many products produced by print finishers.

<u>ADHESIVE PROPERTIES</u>		<u>20 min</u>	<u>24 hour</u>	<u>TEST METHOD</u>
Adhesion to steel	open side	1300 N/M	1300 N/M	BS7116:1990 – D
	closed side	1300 N/M	1300 N/M	BS7116:1990 – D
Shear Adhesion	open side	100 hours		BS7116:1990 – B
	closed side	100 hours		BS7116:1990 – B
Release of liner from adhesive				BS7116:1990 – E

THERMAL PROPERTIES

Recommended service temperature range

0°C to 70°C

Application temperature range

10°C to 50°C

SHELF LIFE

Minimum shelf life when stored in closed cartons away from direct sources of heat and sunlight

2 Years

HEALTH AND SAFETY INFORMATION

Adhesive conforms to F&DA Regulation 175-105

Note: Values given are average and should not be taken as a maximum or minimum for specific purpose. Specific technical details, application information and suggested uses are available. All data relating to suitable uses and description information concerning this product are compiled for research and are believed to be reliable but are provided for guidance only.

For More Details, please contact Tape Range Ltd – www.adhesivetape.net – Tel: 0161 351 9515